CARDIO with HUR SmartTouch Recumbent Ergometer

Easy positioning and entrance for all test subjects





Highlights

Reliable and reproducible stress tests

The experience of professionals who calibrate many ergometers showed that the Lode ergometers are the most reliable across the complete workload and rpm range and still within specifications even after many years of intensive use.

Maximum patient weight 300 kg

Patients of up to 300 kg body weight can be tested on this ergometer

Electric adjustable seat positioning

The seat of the ergometer can be adjusted electrically.

Adjustable back support with flat CPR position

The back support can be adjusted in order to guarantee the most comfort position for each individual test subject. In case of emergency, the back support can be adjusted to a flat CPR position in 1 action.

Suitable for all body postures

This ergometer is designed to offer comfort ergometry testing for test subject of all body postures.





CARDIO with HUR SmartTouch Recumbent Ergometer



Easy positioning and entrance for all test subjects

The latest design of this Recumbent ergometer offers ergometry testing for test subjects up to 300 kg! Easy positioning and access is guaranteed by the electric adjustable seat and adjustable handgrips next to the seat. Moreover, the back support can be adjusted. For safety reasons, the back support can be positioned flat in CPR position in 1 step. This ergometer can be used for elderly and obese people and in rehabilitation settings. It is is standard supplied with a programmable 7" control unit with touch screen, heartrate option and a bottle holder.

A USB A-B cable for service purposes only will be standard delivered with the product.

To connect to the HUR smart touch server you only need a UTP network cable

Features



Extreme low start up load

The extreme low start-up load of 7 watts and the adjustability in small steps of 1 watt **Watt** make this ergometer perfectly suitable for many different applications. The standard control unit shows multiple ergometry parameters and you can determine your specific default setting and start-up menu.



Accurate over a long period of time

The Lode ergometers are supplied with an electro-magnetic braking mechanism of Lanooy (eddy current). The biggest advantage of this braking system compared to a friction braking system is the absolute accuracy and the accuracy over time. Moreover, friction braking systems have more wearing parts.



Small adjustment steps

The workload of the Lode ergometers is adjustable in steps of only 1 watt. Depending **Watt** on your wishes, the test operator or the test subject can adjust the workload. The steps of 1 watt are possible in the manual mode as well as within protocols.



Easy step-through

The ergometer has a comfortable stepthrough: a must for people who are not so mobile!



CARDIO with HUR SmartTouch Recumbent Ergometer

Easy positioning and entrance for all test subjects

CARDIO with HUR SmartTouch Recumbent Ergometer can a.o be extended with the following options:

HR Belt Bluetooth Smart

Smart connection

Pedals with clicker straps

User friendly pedals



Partnumber: 928639



Partnumber: 401111



CARDIO with HUR SmartTouch Recumbent Ergometer



Easy positioning and entrance for all test subjects **Specifications**

Workload			User Interface	
Workload range fixed torque	0,1 - 70 Nm		English user interface	✓
Minimum load	7 W		Chinese user interface	✓
Maximum peak load	1000 W		Croatian user interface	✓
Minimum load increments	1 W		Czech user interface	✓
Maximum continuous load	750 W		Danish user interface	✓
Hyperbolic workload control	~		Dutch user interface	✓
Linear workload control	~		Finnish user interface	✓
Fixed torque workload control	~		French user interface	✓
Maximum rpm independent constant load	150 rpm		German user interface	✓
Minimum rpm independent constant load	30 rpm		Greek user interface	✓
Optional heart rate controlled workload	~		Hungarian user interface	✓
Electromagnetic "eddy current" braking system	~		Italian user interface	✓
Dynamic calibration	~		Japanese user interface	✓
Accuracy			Korean user interface	✓
Workload accuracy from 7 to 100 W	3 W		Latvian user interface	✓
Workload accuracy from 100 to 500 W	3 %		Lithuanian user interface	✓
Workload accuracy from 500 to 1000 W	5 %		Norwegian user interface	✓
Comfort			Polish user interface	✓
Minimum leg length user (incl. adjustable pedals)	650 mm	25.6 inch	Portugese user interface	✓
Allowed user weight	300 kg	661.4 lbs	Romanian user interface	✓
Adjustability range seat	400 mm	15.7 inch	Russian user interface	~
			Spanish user interface	~
			Swedish user interface	~
			Turkish user interface	~
			Ukrainian user interface	~
			Readout Distance	~
			Readout RPM	~
			Readout Heartrate	~
			Readout target HR	~
			Readout Energy	~
			Readout Torque	~
			Readout Time	~
			Readout Power	~
			Set Display	~
			Set Resistance	~
			Set P-Slope	~
			Set Mode	✓
			Manual operation mode	~
			Terminal operation mode	<i>y</i>
				*

Selfdesigned protocol operation mode



CARDIO with HUR SmartTouch Recumbent Ergometer

lla

lla



Easy positioning and entrance for all test subjects

_					
Co	nr	۱۵	ct	11/	ıπ.

Programmable Control Unit with 7" Touchscreen	~	
Dimensions		
Product length (cm)	183 cm	72 inch
Product width (cm)	79 cm	31.1 inch
Product height	149 cm	58.7 inch
Product weight	107 kg	235.9 lbs
Power requirements		
V AC	100 - 240 V	
Phases	1	
Frequency	50/60 Hz	
Power consumption	160 W	
Power cord IEC 60320 C13 with CEE 7/7 plug	~	
Power cord NEMA	×	
Standards & Safety		
IEC 60601-1:2012	~	
ISO 13485:2016 compliant	~	
ISO 9001:2015 compliant	~	

Order info

Certification

Partnumber: 969906

CE class Im according to MDD93/42/EEC
CE class of product with optional SpO2

CE class of product with optional BPM

CB according to IECEE CB - pending



^{*}Specifications are subject to change without notice.